

GKAP1 Antibody



PACO36418

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Golgi apparatus, signal transduction.

Gene ID:

GKAP1

Uniprot

Q5VSY0

Synonyms:

G kinase-anchoring protein 1 (cGMP-dependent protein kinase-anchoring protein of 42 kDa), GKAP1, GKAP42


Immunogen:

Recombinant Human G kinase-anchoring protein 1 protein (21-366AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36418 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO36418 at dilution of 1:100.