

TBCEL Antibody



PACO36446

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Acts as a regulator of tubulin stability.

Gene ID:

TBCEL

Uniprot

Q5QJ74

Synonyms:

Tubulin-specific chaperone cofactor E-like protein (EL) (Leucine-rich repeat-containing protein 35), TBCEL, LRRC35


Immunogen:

Recombinant Human Tubulin-specific chaperone cofactor E-like protein (1-424AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TBCEL antibody at 12 μ g/ml. Lane 1: Rat gonadal tissue. Lane 2: PC-3 whole cell lysate. Lane 3: Rat heart tissue. Lane 4: 293T whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 49 kDa. Observed band size: 49 kDa.