

TBC1D9B Antibody



PACO36470

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May act as a GTPase-activating protein for Rab family protein(s).

Gene ID:

TBC1D9B

Uniprot

Q66K14

Synonyms:

TBC1 domain family member 9B, TBC1D9B, KIAA0676


Immunogen:

Recombinant Human TBC1 domain family member 9B protein (724-1023AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO36470 at dilution of 1:100.



Immunofluorescent analysis of A431 cells using PACO36470 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).