

C6orf141 Antibody




PACO36486


Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: C6orf141
Source: Rabbit	Uniprot Q5SZD1
Isotype: IgG	Synonyms: Uncharacterized protein C6orf141, C6orf141
Applications: ELISA, IHC	Immunogen: Recombinant Human Uncharacterized protein C6orf141 protein (1-237AA).
Recommended dilutions: ELISA:1:2000-1:10000, IHC:1:20-1:200	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO36486 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO36486 at dilution of 1:100.