

NXNL2 Antibody, Biotin conjugated



PACO36501

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the maintenance of both the function and the viability of sensory neurons, including photoreceptors and olfactory neurons.

Gene ID:

NXNL2

Uniprot

Q5VZ03

Synonyms:

Nucleoredoxin-like protein 2 (Rod-derived cone viability factor 2) (RdCVF2), NXNL2, C9orf121

Immunogen:

Recombinant Human Nucleoredoxin-like protein 2 protein (1-135AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A