

ATL3 Antibody, FITC conjugated



PACO36524

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GTPase tethering membranes through formation of trans-homooligomers and mediating homotypic fusion of endoplasmic reticulum membranes. Functions in endoplasmic reticulum tubular network biogenesis.

Gene ID:

ATL3

Uniprot

Q6DD88

Synonyms:

Atlastin-3 (EC 3.6.5. -), ATL3

Immunogen:

Recombinant Human Atlastin-3 protein (1-187AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A