

# ATL3 Antibody, Biotin conjugated



PACO36525

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GTPase tethering membranes through formation of trans-homooligomers and mediating homotypic fusion of endoplasmic reticulum membranes. Functions in endoplasmic reticulum tubular network biogenesis.

**Gene ID:**

ATL3

**Uniprot**

Q6DD88

**Synonyms:**

Atlastin-3 (EC 3.6.5. -), ATL3

**Immunogen:**

Recombinant Human Atlastin-3 protein (1-187AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A