

LGALS9C Antibody




PACO36526

Product Information

Size:	Protein Background:
50ug	Binds galactosides.
Reactivity:	Gene ID:
Human, Mouse	LGALS9C
Source:	Uniprot
Rabbit	Q6DKI2
Isotype:	Synonyms:
IgG	Galectin-9C (Gal-9C) (Galectin-9-like protein B), LGALS9C
Applications:	Immunogen:
ELISA, WB	Recombinant Human Galectin-9C protein (1-356AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:1000-1:5000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: LGALS9C antibody at 4 μ g/ml. Lane 1: Mouse stomach tissue. Lane 2: Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 40 kDa. Observed band size: 22 kDa.