

LDLRAD1 Antibody



PACO36542

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

LDLRAD1

Uniprot

Q5T700

Synonyms:

Low-density lipoprotein receptor class A domain-containing protein 1, LDLRAD1


Immunogen:

Recombinant Human Low-density lipoprotein receptor class A domain-containing protein 1 protein (71-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO36542 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).