

C9orf57 Antibody



PACO36562

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C9orf57

Uniprot

Q5W0N0

Synonyms:

Uncharacterized protein C9orf57, C9orf57


Immunogen:

Recombinant Human Uncharacterized protein C9orf57 protein (51-127AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO36562 at dilution of 1:100.