

PACO36566

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transport of lipophilic molecules, possible pheromone-carrier.

**Gene ID:**

Lcn3

**Uniprot**

Q62471

**Synonyms:**

Vomer nasal secretory protein 1 (Lipocalin-3) (Vomer nasal secretory protein I) (VNSP I), Lcn3

**Immunogen:**

Recombinant Mouse Vomer nasal secretory protein 1 protein (19-182AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A