

ZKSCAN2 Antibody



PACO36570

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZKSCAN2

Uniprot

Q63HK3

Synonyms:

Zinc finger protein with KRAB and SCAN domains 2 (Zinc finger protein 694), ZKSCAN2, ZNF694


Immunogen:

Recombinant Human Zinc finger protein with KRAB and SCAN domains 2 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO36570 at dilution of 1:100.