

# Cyp3a13 Antibody



PACO36578

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Can activate aflatoxin B1 to a genotoxic product.

**Gene ID:**

Cyp3a13

**Uniprot**

Q64464

**Synonyms:**

Cytochrome P450 3A13 (EC 1.14.14.1) (CYP11A13), Cyp3a13, Cyp3a-13

**Immunogen:**

Recombinant Mouse Cytochrome P450 3A13 protein (1-503AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A