

# FAM69A Antibody



PACO36594


---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	FAM69A
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q5T7M9
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Protein FAM69A, FAM69A
<b>Applications:</b>	<b>Immunogen:</b>
ELISA, WB	Recombinant Human Protein FAM69A protein (49-162AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
ELISA:1:2000-1:10000, WB:1:1000-1:5000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western blot. All lanes: FAM69A antibody at 8 $\mu$ g/ml + HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 50, 19 kDa. Observed band size: 50 kDa.