

RPL7L1 Antibody



PACO36638

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytosolic large ribosomal subunit, poly(A) RNA binding, structural constituent of ribosome, cytoplasmic translation, maturation of LSU-rRNA from tricistronic rRNA transcript (SSU-rRNA, 5.8S rRNA, LSU-rRNA).

Gene ID:

RPL7L1

Uniprot

Q6DKI1

Synonyms:

60S ribosomal protein L7-like 1 (Large ribosomal subunit protein uL30-like 1), RPL7L1


Immunogen:

Recombinant Human 60S ribosomal protein L7-like 1 protein (1-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO36638 at dilution of 1:100.