

# LYPLAL1 Antibody, HRP conjugated



PACO36659

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Has depalmitoylating activity toward KCNMA1. Does not exhibit phospholipase nor triacylglycerol lipase activity, able to hydrolyze only short chain substrates due to its shallow active site.

**Gene ID:**

LYPLAL1

**Uniprot**

Q5VWZ2

**Synonyms:**

Lysophospholipase-like protein 1 (EC 3.1.2. -), LYPLAL1

**Immunogen:**

Recombinant Human Lysophospholipase-like protein 1 protein (2-237AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A