

Arachin Ahy-3 Antibody



PACO36666

Product Information

Size:

50ug

Reactivity:

Peanut

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Q647H2

Synonyms:

Arachin Ahy-3 [Cleaved into: Arachin Ahy-3 chain alpha; Arachin Ahy-3 chain beta]

Immunogen:

Recombinant Peanut Arachin Ahy-3 protein (299-478AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A