

PACO36686

---

## Product Information

**Size:**

50ug

**Reactivity:**

Staphylococcus aureus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transfers the 4'-phosphopantetheine moiety from coenzyme A to a Ser of acyl-carrier-protein.

**Gene ID:**

acpS

**Uniprot**

Q6GF02

**Synonyms:**

Holo-[acyl-carrier-protein] synthase (Holo-ACP synthase) (EC 2.7.8.7) (4'-phosphopantetheinyl transferase AcpS), acpS, dpj

**Immunogen:**

Recombinant Staphylococcus aureus Holo-[acyl-carrier-protein] synthase protein (1-119AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A