

ZNF767P Antibody




PACO36730


Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	ZNF767P
Source:	Uniprot
Rabbit	Q75MW2
Isotype:	Synonyms:
IgG	Protein ZNF767 (Zinc finger protein 767 pseudogene), ZNF767P, ZNF767
Applications:	Immunogen:
ELISA, IHC	Recombinant Human Protein ZNF767 protein (1-118AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO36730 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO36730 at dilution of 1:100.