

TANGO2 Antibody



PACO36742

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TANGO2

Uniprot

Q6ICL3

Synonyms:

Transport and Golgi organization protein 2 homolog, TANGO2, C22orf25

Immunogen:

Recombinant Human Transport and Golgi organization protein 2 homolog protein (1-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A