

SLC30A6 Antibody, FITC conjugated



PACO36756

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Zinc-efflux transporter which allocates the cytoplasmic zinc to the trans-Golgi network (TGN) as well as the vesicular compartment.

Gene ID:

SLC30A6

Uniprot

Q6NXT4

Synonyms:

Zinc transporter 6 (ZnT-6) (Solute carrier family 30 member 6), SLC30A6, ZNT6

Immunogen:

Recombinant Human Zinc transporter 6 protein (246-461AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A