

AGTRAP Antibody, FITC conjugated



PACO36764

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Appears to be a negative regulator of type-1 angiotensin II receptor-mediated signaling by regulating receptor internalisation as well as mechanism of receptor desensitization such as phosphorylation. Induces also a decrease in cell proliferation and angiotensin II-stimulated transcriptional activity.

Gene ID:

AGTRAP

Uniprot

Q6RW13

Synonyms:

Type-1 angiotensin II receptor-associated protein (AT1 receptor-associated protein), AGTRAP, ATRAP

Immunogen:

Recombinant Human Type-1 angiotensin II receptor-associated protein (108-152AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A