

PQLC2 Antibody



PACO36786

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Amino acid, transporter that specifically mediates the pH-dependent export of the cationic amino acid, arginine, histidine and lysine from lysosomes.

Gene ID:

PQLC2

Uniprot

Q6ZP29

Synonyms:

Lysosomal amino acid, transporter 1 homolog (PQ-loop repeat-containing protein 2), PQLC2


Immunogen:

Recombinant Human Lysosomal amino acid, transporter 1 homolog protein (91-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO36786 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).