

KCTD1 Antibody



PACO36794

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May repress the transcriptional activity of AP-2 family members, including TFAP2A, TFAP2B and TFAP2C to various extent.

Gene ID:

KCTD1

Uniprot

Q719H9

Synonyms:

BTB/POZ domain-containing protein KCTD1 (Potassium channel tetramerization domain-containing protein 1), KCTD1, C18orf5

Immunogen:

Recombinant Human BTB/POZ domain-containing protein KCTD1 protein (1-257AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A