

MORN3 Antibody




PACO36838

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse	Gene ID: MORN3
Source: Rabbit	Uniprot Q6PF18
Isotype: IgG	Synonyms: MORN repeat-containing protein 3, MORN3
Applications: ELISA, WB	Immunogen: Recombinant Human MORN repeat-containing protein 3 protein (1-240AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: MORN3 antibody at 2 μ g/ml + Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 28, 19 kDa. Observed band size: 28 kDa.