

# TTC33 Antibody



PACO36842

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:****Gene ID:**

TTC33

**Uniprot**

Q6PID6

**Synonyms:**

Tetratricopeptide repeat protein 33 (TPR repeat protein 33) (Osmosis-responsive factor), TTC33

**Immunogen:**


Recombinant Human Tetratricopeptide repeat protein 33 protein (1-262AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36842 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO36842 at dilution of 1:100.