

# LRRN1 Antibody, FITC conjugated



PACO36856

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LRRN1

**Uniprot**

Q6UXK5

**Synonyms:**

Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1), LRRN1, KIAA1497

**Immunogen:**

Recombinant Human Leucine-rich repeat neuronal protein 1 protein (26-631AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A