

LRRN1 Antibody, Biotin conjugated



PACO36857

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LRRN1

Uniprot

Q6UXK5

Synonyms:

Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1), LRRN1, KIAA1497

Immunogen:

Recombinant Human Leucine-rich repeat neuronal protein 1 protein (26-631AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A