

GRHL2 Antibody, Biotin conjugated



PACO36881

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GRHL2

Uniprot

Q6ISB3

Synonyms:

Grainyhead-like protein 2 homolog (Brother of mammalian grainyhead) (Transcription factor CP2-like 3), GRHL2, BOM TFCP2L3

Immunogen:

Recombinant Human Grainyhead-like protein 2 homolog protein (1-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A