

FAM71F2 Antibody




PACO36886


Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	FAM71F2
Source:	Uniprot
Rabbit	Q6NXP2
Isotype:	Synonyms:
IgG	Protein FAM71F2 (Protein FAM137B), FAM71F2, FAM137B
Applications:	Immunogen:
ELISA, IHC	Recombinant Human Protein FAM71F2 protein (1-300AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO36886 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36886 at dilution of 1:100.