

LYPD4 Antibody



PACO36902

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LYPD4

Uniprot

Q6UWN0

Synonyms:

Ly6/PLAUR domain-containing protein 4, LYPD4

Immunogen:

Recombinant Human Ly6/PLAUR domain-containing protein 4 protein (46-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A