

GBP6 Antibody



PACO36914

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds GTP, GDP and GMP.

Gene ID:

GBP6

Uniprot

Q6ZLN6

Synonyms:

Guanylate-binding protein 6 (GTP-binding protein 6) (GBP-6) (Guanine nucleotide-binding protein 6), GBP6


Immunogen:

Recombinant Human Guanylate-binding protein 6 protein (401-600AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO36914 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36914 at dilution of 1:100.