

RUFY4 Antibody



PACO36918

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

RUFY4

Uniprot

Q6ZNE9

Synonyms:

RUN and FYVE domain-containing protein 4, RUFY4


Immunogen:

Recombinant Human RUN and FYVE domain-containing protein 4 protein (174-473AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36918 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung tissue using PACO36918 at dilution of 1:100.