

Polygalacturonase Antibody



PACO36930

Product Information

Size:

50ug

Reactivity:

Chamaecyparis obtusa

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Random hydrolysis of (1->4)-alpha-D-galactosiduronic linkages in pectate and other galacturonans.

Gene ID:**Uniprot**

Q7M1E7

Synonyms:

Polygalacturonase (PG) (EC 3.2.1.15) (Major pollen allergen Cha o 2) (Pectinase) (allergen Cha o 2)

Immunogen:

Recombinant Chamaecyparis obtusa Polygalacturonase protein (51-514AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A