

PACO36934

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cell surface-associated protein implicated in virulence by promoting bacterial attachment to both alpha- and beta-chains of human fibrinogen and inducing the formation of bacterial clumps.

Gene ID:

clfB

Uniprot

Q6GDH2

Synonyms:

Clumping factor B (Fibrinogen receptor B) (Fibrinogen-binding protein B), clfB

Immunogen:

Recombinant Staphylococcus aureus Clumping factor B protein (45-542AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A