

# FAM109B Antibody



PACO36942

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a role in endocytic trafficking. Required for receptor recycling from endosomes, both to the trans-Golgi network and the plasma membrane.

**Gene ID:**

FAM109B

**Uniprot**

Q6ICB4

**Synonyms:**

Sesquipedalian-2 (Ses2) (27 kDa inositol polyphosphate phosphatase interacting protein B) (IPIP27B), FAM109B

**Immunogen:**

Recombinant Human Sesquipedalian-2 protein (1-259AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A