

ETFRF1 Antibody



PACO36946

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as a regulator of the electron transfer flavoprotein by promoting the removal of flavin from the ETF holoenzyme (composed of ETFA and ETFB).

Gene ID:

ETFRF1

Uniprot

Q6IPR1

Synonyms:

Electron transfer flavoprotein regulatory factor 1 (LYR motif-containing protein 5), ETFRF1, LYRM5


Immunogen:

Recombinant Human Electron transfer flavoprotein regulatory factor 1 protein (1-88AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO36946 at dilution of 1:100.