

OGFOD2 Antibody



PACO36954

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse	Gene ID: OGFOD2
Source: Rabbit	Uniprot Q6N063
Isotype: IgG	Synonyms: 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 2 (EC 1.14.11. -), OGFOD2
Applications: ELISA, WB	Immunogen: Recombinant Human 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 2 protein (1-186AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: OGFOD2 antibody at 3 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 39, 33, 15, 21 kDa. Observed band size: 39 kDa.