

# NANOGP8 Antibody



PACO36958

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

May act as a transcription regulator. When overexpressed, promotes entry of cells into S phase and cell proliferation.

**Gene ID:**

NANOGP8

**Uniprot**

Q6NSW7

**Synonyms:**

Homeobox protein NANOGP8, NANOGP8

**Immunogen:**


Recombinant Human Homeobox protein NANOGP8 protein (1-305AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Rat liver tissue, Mouse liver tissue. All lanes: NANOGP8 antibody at 2 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 35 kDa. Observed band size: 35 kDa.