

# FAM117B Antibody



PACO36966

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

FAM117B

**Uniprot**

Q6P1L5

**Synonyms:**

Protein FAM117B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein), FAM117B, ALS2CR13

**Immunogen:**


Recombinant Human Protein FAM117B protein (1-200AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human glioma using PACO36966 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO36966 at dilution of 1:100.