

TUBA3E Antibody, FITC conjugated



PACO36972

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Gene ID:

TUBA3E

Uniprot

Q6PEY2

Synonyms:

Tubulin alpha-3E chain (Alpha-tubulin 3E) [Cleaved into: Detyrosinated tubulin alpha-3E chain], TUBA3E

Immunogen:

Recombinant Human Tubulin alpha-3E chain protein (1-250AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A