

PACO36974

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Potently stimulates the neddylation of cullin components of SCF-type E3 ubiquitin ligase complexes from the NEDD8-conjugating E2 enzyme UBC12. Neddylation of cullins play an essential role in the regulation of SCF-type complexes activity.

Gene ID:

DCUN1D2

Uniprot

Q6PH85

Synonyms:

DCN1-like protein 2 (DCUN1 domain-containing protein 2) (Defective in cullin neddylation protein 1-like protein 2), DCUN1D2, C13orf17 DCUN1L2

Immunogen:

Recombinant Human DCN1-like protein 2 protein (1-259AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A