

KIRREL2 Antibody, HRP conjugated



PACO36983

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May regulate basal insulin secretion.

Gene ID:

KIRREL2

Uniprot

Q6UWL6

Synonyms:

Kin of IRRE-like protein 2 (Kin of irregular chiasm-like protein 2) (Nephrin-like protein 3), KIRREL2, NEPH3

Immunogen:

Recombinant Human Kin of IRRE-like protein 2 protein (21-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A