

# DENND5A Antibody, HRP conjugated



PACO37003

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Guanine nucleotide exchange factor (GEF) which may activate RAB6A and RAB39A and/or RAB39B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form.

**Gene ID:**

DENND5A

**Uniprot**

Q6IQ26

**Synonyms:**

DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1), DENND5A, KIAA1091 RAB6IP1

**Immunogen:**

Recombinant Human DENN domain-containing protein 5A protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A