

FOXR2 Antibody, Biotin conjugated



PACO37025

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleus, RNA polymerase II transcription factor activity, sequence-specific DNA binding, anatomical structure morphogenesis, cell differentiation, regulation of transcription from RNA polymerase II promoter.

Gene ID:

FOXR2

Uniprot

Q6PJQ5

Synonyms:

Forkhead box protein R2 (Forkhead box protein N6), FOXR2, FOXN6

Immunogen:

Recombinant Human Forkhead box protein R2 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A