

CFAP65 Antibody, HRP conjugated



PACO37035

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in sperm motility.

Gene ID:

CFAP65

Uniprot

Q6ZU64

Synonyms:

Cilia- and flagella-associated protein 65 (Coiled-coil domain-containing protein 108), CFAP65, CCDC108

Immunogen:

Recombinant Human Cilia- and flagella-associated protein 65 protein (2-164AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A