

# FKBP9P1 Antibody



PACO37042

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FKBP9P1

**Uniprot**

Q75LS8

**Synonyms:**

Putative FK506-binding protein 9-like protein (FK506-binding protein 9-like protein pseudogene), FKBP9P1, FKBP9L

**Immunogen:**

Recombinant Human Putative FK506-binding protein 9-like protein (1-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A