

DHRS7B Antibody



PACO37054

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Putative oxidoreductase.

Gene ID:

DHRS7B

Uniprot

Q6IAN0

Synonyms:

Dehydrogenase/reductase SDR family member 7B (EC 1.1. -. -) (Short-chain dehydrogenase/reductase family 32C member 1), DHRS7B, SDR32C1

Immunogen:

Recombinant Human Dehydrogenase/reductase SDR family member 7B protein (39-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A