

# ASPHD2 Antibody, FITC conjugated



PACO37060

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as 2-oxoglutarate-dependent dioxygenase.

**Gene ID:**

ASPHD2

**Uniprot**

Q6ICH7

**Synonyms:**

Aspartate beta-hydroxylase domain-containing protein 2 (EC 1.14.11. -), ASPHD2

**Immunogen:**

Recombinant Human Aspartate beta-hydroxylase domain-containing protein 2 protein (80-369AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A