

PAXIP1 Antibody, FITC conjugated



PACO37120

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PAXIP1

Uniprot

Q6ZW49

Synonyms:

PAX-interacting protein 1 (PAX transactivation activation domain-interacting protein), PAXIP1, PAXIP1L PTIP

Immunogen:

Recombinant Human PAX-interacting protein 1 protein (613-812AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A